








A PATIENT OVERVIEW

Your treatment options for
Stress Urinary Incontinence (SUI)

Overview of treatment options for male stress urinary incontinence (SUI)

Individual results may vary and are not predictive of results in other cases. Consult with your doctor to decide the option most suited for you.

	Solution satisfaction and attributes	Some possible side effects	Characteristics	Typical duration of use
 <p>Male sling</p> <p>A sling made of soft mesh implanted inside the body to reposition the urethra and provide support to surrounding sphincter muscle. This can help to keep the urethra closed, especially when coughing, sneezing, and lifting.</p>	<ul style="list-style-type: none"> • 79.5% surveyed would have the AdvVance™ XP Male Sling procedure again³ • 88% surveyed would recommend the AdvVance XP Male Sling to a friend³ • More than 10 years of clinical use • Assists natural continence function • Surgical treatment 	<ul style="list-style-type: none"> • Device failure • Urinary retention • Postoperative pain • Irritation at the wound site • Foreign body response 	<ul style="list-style-type: none"> • Long-term treatment option • Sling acts as a "hammock," repositioning and supporting the urethra to help restore normal bladder control • Made of synthetic mesh • Placed entirely inside the body, making it undetectable to others • No mechanical components that require operation 	<ul style="list-style-type: none"> • Permanent implant • Demonstrated cured and improved rates ranging from 63% to 83% at 5 years in patients with mild to moderate incontinence^{4,5} • All alternate options still available if the sling loses effectiveness
 <p>Artificial urinary sphincter</p> <p>The artificial urinary sphincter (AUS) is placed inside the body. A fluid-filled cuff keeps the urethra closed and a pump in the scrotum allows urination on demand. The AUS is designed for all levels of SUI following prostate procedures and is considered the gold standard treatment.⁶⁻⁸</p>	<ul style="list-style-type: none"> • 94.6% of patients* are satisfied with their AMS 800™ Artificial Urinary Sphincter (AUS)⁹ • 94% or more would recommend the AMS 800 AUS procedure to a friend¹⁰ • More than 40 years of clinical use of the AUS, with over 200,000 systems sold worldwide • Placed entirely inside the body, it is undetectable to others • Natural-acting elimination • Surgical solution • Functions all day and night <p>*Non-Urethroplasty patients</p>	<ul style="list-style-type: none"> • Device malfunction or failure, which may require revision surgery • Erosion of the urethra in the cuff area • Urinary retention • Postoperative pain 	<ul style="list-style-type: none"> • Long-term treatment option • Three-part urinary control system contained completely in the body • The control pump is implanted in the scrotum • The inflatable cuff is placed around the urethra • The fluid-filled pressure regulating balloon (PRB) is usually implanted in a natural, open space next to the bladder • Designed to treat male SUI following prostate procedures • Mimics a healthy sphincter, allowing patients to urinate when desired • Requires good cognitive ability and manual dexterity • AUS is considered a gold standard treatment for male SUI⁶⁻⁸ 	<ul style="list-style-type: none"> • Ongoing, with device replacement only as needed • 90% of patients are satisfied with their AMS 800 AUS device long-term^{11*} <p>**10+ years</p>
 <p>Absorbent products</p> <p>Incontinence pads are commonly used to help absorb urine leakage, protect the skin and block odor. Other absorbent products include diapers and undergarments. These products can offer odor control; some are disposed after a single use; others may be cleaned and reused.</p>	<ul style="list-style-type: none"> • Appropriate for everyday activities • Comfortable when dry • Good for overnight • Easy to use 	<ul style="list-style-type: none"> • Requires changing, up to several times per day • Can leak • Uncomfortable when wet • Expensive over time 	<ul style="list-style-type: none"> • Pads, liners, or undergarments typically have a waterproof backing • Some products contain a breathable plastic film that helps reduce skin irritation 	<ul style="list-style-type: none"> • Requires frequent changing on an hourly/daily basis
 <p>External collection devices</p> <p>Condom catheters or a penis pouch are placed on the penis so urine can flow into a drainage bag. Urine collection bags are strapped to the body underneath clothing during the day and may rest bedside at night.</p>	<ul style="list-style-type: none"> • Discrete option for long periods of time • Keeps skin dry • Avoids odor issues • Convenient for storage and travel 	<ul style="list-style-type: none"> • Device can move or loosen resulting in leaks • Irritative or allergic reactions can occur based on material¹² • Long-term use may result in urinary tract infections¹² 	<ul style="list-style-type: none"> • Typically made of silicone or latex • Uses adhesive to stay in place 	<ul style="list-style-type: none"> • Designed to be worn 24/7 and changed daily
 <p>Penile clamp</p> <p>A penile clamp uses a hinged, rigid frame that supports two pads and a locking mechanism. It controls leakage by applying constant pressure upon the penis.</p>	<ul style="list-style-type: none"> • Can stop leaking • Simple to remove • Can be washed and reused 	<ul style="list-style-type: none"> • Can only be worn for short periods of time¹³ • May be described as uncomfortable and painful • Swelling¹³ • Strictures of the urethra¹³ • Development of skin breakdown¹³ 	<ul style="list-style-type: none"> • Most have a padded, flexible soft interior of foam that conforms to fit the penis • Exterior of the clamp may be metal or plastic 	<ul style="list-style-type: none"> • MUST be removed every 1-2 hours to allow for urination¹³ • Should not be worn 24 hours, 7 days a week¹³

Treatment options for male stress urinary incontinence (SUI)

Urinary incontinence is a common complication following prostate cancer surgery, with 15 % to 90 % of men experiencing incontinence at one year after surgery.¹ Incontinence, also known as bladder leakage, can stop you from living the active life you once lived and make you feel alone and isolated. For many men, leaking urine may lead to feelings of embarrassment and depression.²

With a thorough evaluation of your bladder leakage symptoms, your doctor may recommend changes to your daily routine, such as avoiding caffeine and alcohol. You may modify your fluid intake and follow a voiding schedule to try to train your bladder. Your doctor may also recommend pelvic floor exercises to tighten the sphincter muscles.

If bladder leakage persists, some men cope by using absorbent pads and protective undergarments. Others use special devices, such as catheters or penile clamps, to try to prevent urine leakage.

Long-term treatments for bladder leakage include the AdVance™ XP Male Sling System and the AMS 800™ Artificial Urinary Sphincter (AUS). Both options are placed entirely within the body for a discrete solution for SUI. If bladder leakage becomes a burden and affects your happiness, you can talk to a urologist specializing in male stress incontinence about options to restore bladder control, and find the right treatment for you.

Take the next step

Don't let stress urinary incontinence (SUI) stop you from living your best life. Learn about SUI treatment options, find a specialist, and hear from patients who found effective SUI solutions.

Important Risk Information

AdVance™ XP Male Sling System: Your doctor is your best source for information on the risks and benefits of the AdVance™ XP Male Sling System. Talk to your doctor for a complete listing of risks, warnings and important safety information.

The AdVance™ XP Male Sling System is intended for the treatment of male stress urinary incontinence (SUI). Potential risks may include inability to urinate (urinary retention), return to incontinence, infection, erosion, device migration, pelvic organ dysfunction, bleeding, and pain. These devices are contraindicated for patients with urinary tract infections, blood coagulation disorders, a compromised immune system or any other condition that would compromise healing, with renal insufficiency, and upper urinary tract relative obstruction.

AMS 800™ Artificial Urinary Sphincter: Your doctor is your best source for information on the risks and benefits of the AMS 800™ Artificial Urinary Sphincter. Talk to your doctor for a complete listing of risks, warnings and important safety information.

The AMS 800™ Artificial Urinary Sphincter is intended for use in the treatment of male stress urinary incontinence (intrinsic sphincter deficiency) following prostate surgery.

Men with diabetes, spinal cord injuries or skin infections may have an increased risk of infection. Some AMS 800 devices contain an antibiotic (InhibiZone™ Antibiotic Surface Treatment). The device may not be suited for patients who are allergic to the antibiotics contained within the device (rifampin, minocycline or other tetracyclines) or have systemic lupus.

Potential risks may include: device malfunction/failure leading to additional surgery, wearing away/loss of tissue (device/tissue erosion), inability to urinate (urinary retention), infection, and pain/soreness.

References

1. Gupta S, Ding L, Granieri M, Le N-B, Peterson AC. Utilization of surgical procedures and racial disparity in the treatment of urinary incontinence after prostatectomy. *Neurourol Urodyn*. 2016;35:733–737.
2. Wu S, Wu F. Association of urinary incontinence with depression among men: a cross-sectional study. *BMC Public Health*. 2023;23:944.
3. Hüscher T, Kretschmer A, Thomsen F, et al. The AdVance and AdVanceXP male sling in urinary incontinence: is there a difference? *World J Urol*. 2018;36:1657–1662.
4. Mumm JN, Klehr B, Rodler S, et al. Five-year results of a prospective multicenter trial: AdVance XP for postprostatectomy-incontinence in patients with favorable prognostic factors. *Urol Int*. 2021;105:421–427.
5. Rizvi IG, Ravindra P, Pipe M, Sohawon R, King T, Belal M. The AdVance™ male sling: does it stand the test of time? *Scand J Urol*. 2021;55:155–160.
6. Biardeau X, Aharoni S; AUS Consensus Group, et al. Artificial Urinary Sphincter: Report of the 2015 Consensus Conference. *Neurourol Urodyn*. 2016;35(Suppl 2):S8–S24.
7. Montague DK. Artificial urinary sphincter: long-term results and patient satisfaction. *Adv Urol*. 2012;2012:835290.
8. James MH, McCammon KA. Artificial urinary sphincter for post-prostatectomy incontinence: a review. *Int J Urol*. 2014;21:536–543.
9. Sayedahmed K, Ollanas R, Kaftan B, et al. Impact of previous urethroplasty on the outcome after artificial urinary sphincter implantation: a prospective evaluation. *World J Urol*. 2020;38:183–191.
10. Linder BJ, Rivera ME, Ziegelmann MJ, et al. Long-term outcomes following artificial urinary sphincter placement: an analysis of 1082 cases at Mayo Clinic. *Urology*. 2015;86:602–607.
11. Léon P, Chartier-Kastler E, Roupêt M, et al. Long-term functional outcomes after artificial urinary sphincter implantation in men with stress urinary incontinence. *BJU Int*. 2015;115:951–957.
12. Evidence-based guidelines for best practice in urological health care. Male external catheters in adults. Urinary catheter management. European Association of Urology Nurses. March 2016.
13. Penile Clamps/Cuffs for Male Urinary Incontinence. The Simon Foundation for Continence. <https://simonfoundation.org/penile-clamp-urinary-incontinence>. Accessed May 1, 2023.
14. Penile Clamps Q & A. Seek Wellness. www.seekwellness.com/male_incontinence/penile-clampsqa.htm. Accessed July 13, 2017.

This material is for informational purposes only and not meant for medical diagnosis. This information does not constitute medical or legal advice, and Boston Scientific makes no representation regarding the medical benefits included in this information. Boston Scientific strongly recommends that you consult with your physician on all matters pertaining to your health.

All images are the property of Boston Scientific.

All trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician.

Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2024 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-2002909-AA

CE 2797